

Reverse Osmosis Anti-scaling

Product	Chemical nature	Dose
---------	-----------------	------

Inhibitor of scaling for osmosis

INCUS®-RO/3	Product based on acrylic polymers and phosphonic derivatives.	0,5-3 ppm over the flow rate to be treated
INCUS®-RO/40	Product based on acrylic polymers and phosphonic derivatives.	0,5-3 ppm over the flow rate to be treated
INCUS®-RO/50	Product based on acrylic polymers.	0,5-3 ppm over the flow rate to be treated
INCUS®-RO/89	Product based on phosphonic derivatives.	2-6 ppm over the flow rate to be treated.
INCUS®-RO/100	Product based on phosphonic derivatives.	2-6 ppm over the flow rate to be treated.
RESTIN®-RO/75	Based on maleic copolymer.	0,5-3 ppm over the flow rate to be treated
RESTIN®-RO/130	Based on inorganic polyphosphates.	1-10 ppm over the flow rate to be treated.